

FLUXANA

Your supplier for
XRF Application Solutions

HD ELEKTRONIK

Fusion machines

Liquid Analysis Technology



FLUXANA GmbH & Co. KG
Borschelstr. 3, 47551 Bedburg-Hau, Germany
Tel.: +49 (0) 2821 997 32-0
Fax: +49 (0) 2821 997 32 29
E-mail: info@fluxana.de
Web: www.fluxana.com

Amtsgericht Kleve: HR-A 2935, HR-B 8211
Ust-IdNr.: DE 814692564, Steuer-Nr. 116/5755/0442
Finanzamt Kleve

In cooperation with:

HD Elektronik und Elektrotechnik GmbH
Sommerdeich 22, 47533 Kleve, Germany
Tel.: +49 (0) 2821 148 10
Fax: +49 (0) 2821 148 09
E-mail: hde@hdelektronik.de
Web: www.hdelektronik.de

Amtsgericht Kleve: HR-B 2162
Ust-IdNr.: DE 812941185, Steuer-Nr. 116/5707/1869
Finanzamt Kleve



Official agent



NEW! Cups with film
Already prepared


Please ask for a quote.




Thin Films

The preferences of our thin films:

- 25% more pre-cut films per package than other brands
- dust proof anti static packaging
- lot number tracking system
- rolls, pre-cut circles & sheets



Item No.	Description	Thickness μm	Packaging
TF-025	Microporous		Roll 6,4cm x 5,1m
TF-112	Mylar®	12	Roll 7,6cm x 91,4m
TF-115	Ultrapolyester	1,5	Roll 7,6cm x 91,4m
TF-125	Mylar®	2,5	Roll 7,6cm x 91,4m
TF-125-345	Mylar®	2,5	Square 7,6cm x 7,6cm (500 pcs.)
TF-12560	PolyRep® (replacement for polycarbonate)	Mylar® 2,5 + PP 6,0	2 Rolls a 7,6cm x 91,4m
TF-1335	Hostaphan®	3,5	Roll 7,6cm x 91,4m
TF-1335-345	Hostaphan®	3,5	Square 7,6cm x 7,6cm (500 pcs.)
TF-135	Mylar®	3,5	Roll 7,6cm x 91,4m
TF-135-255	Mylar®	3,5	Circle 6,4cm (500 pcs.)
TF-135-345	Mylar®	6,0	Square 7,6cm x 7,6cm (500 pcs.)
TF-160	Mylar®	6,0	Roll 7,6cm x 91,4m
TF-160-255	Mylar®	6,0	Circle 6,4cm (500 pcs.)
TF-160-345	Mylar®	6,0	Square 7,6cm x 7,6cm (500 pcs.)
TF-212	Polypropylene	12,0	Roll 7,6cm x 91,4m
TF-212-345	Polypropylene	12,0	Square 7,6cm x 7,6cm (500 pcs.)
TF-240	Polypropylene	4,0	Roll 7,6cm x 91,4m
TF-240-255	Polypropylene	4,0	Circle 6,4cm (500 pcs.)
TF-240-30	Polypropylene	4,0	Circle 7,6cm (100 pcs.)
TF-240-345	Polypropylene	4,0	Square 7,6cm x 7,6cm (500 pcs.)
TF-240-355	Polypropylene	4,0	Circle 8,9cm (500pcs.) (best for SC-3345)
TF-250	Polypropylene	5,0	Roll 7,6cm x 91,4m
TF-260	Polypropylene	6,0	Roll 7,6cm x 91,4m
TF-260-255	Polypropylene	6,0	Circle 6,4cm (500 pcs.)
TF-260-345	Polypropylene	6,0	Square 7,6cm x 7,6cm (500 pcs.)
TF-350	Polycarbonate	5,0	Roll 7,6cm x 91,4m
TF-475	Kapton®	7,5	Roll 7,6cm x 91,4m
TF-475-345	Kapton®	7,5	Square 7,6cm x 7,6cm (500 pcs.)
TF-500	XRF Tape	-	Roll 25mm x 66m
TF-LS-240-255	Polypropylene Low Sulfur	4,0	Circle 6,4cm (500 pcs.)
WZ-0001	Preparation tool for our sample cups	-	Measurement: 56 x 20mm



Professional Cups

The sample cups of FLUXANA are designed for use in XRF. All cups have been optimized and approved by our customers world wide.



We offer different XRF-sample cups to you.

Item No.	Description	OD mm	ID mm	H mm	Vol ml
SC-3332	Cup, ring, cap for: Spectro, XEPOS, XLAB 2000®	31	24	23	9
SC-3335	Cup, ring, cap for Philips/Panalytical®	35	28	35	17
SC-3340	Inner ring, outer ring, cap for: Bruker®, SRS, S4, S8, ARL Optimix®	41	34	25	17
SC-3340T	Inner ring, outer ring, cap for: Bruker® SRS, S4, S8	41	32	36	22
SC-3345	Inner ring, outer ring, cap universal for 45mm	45	37	37	33
SC-3351	Cup, cap universal for 51,5mm holder	51	44	41	48
SC-4031	Cup, ring, single open ended, universal open for 32mm holder	31	24	24	9
SC-4047	Cup, ring, single open ended, for Horiba®	46	35	18	12
SC-4131	Cup, ring, collar single open ended, universal for 32mm holder	32	25	24	9
SC-4140	Cup, ring, collar single open ended, universal for 40mm holder	40	32	23	16
SC-4231	Cup, 2 rings, double open ended, universal for 32mm holder	31	24	24	10
SC-4240	Cup, 2 rings, double open ended, universal for 40mm holder	39	32	23	18
SC-4331	Cup, ring, cap universal for 32mm holder	31	24	26	9
SC-4335	Cup, ring, cap for Philips/Panalytical®	34	27	36	17
SC-4340	Cup, ring, cap universal for 40mm holder	39	31	25	15
SC-4345	Cup, cap for Philips/Panalytical®	45	38	41	35
SC-6232	Cup, 2 rings, double open ended for Spectro XEPOS, XLAB 2000®	32	26	21	11
SC-6332	Cup, ring, cap for Spectro XEPOS, XLAB 2000®	32	26	23	11
SC-7332	Inner cup and cap, designed for Oxford®	34	28	38	16
SC-7332inner	Only inner cup, designed for Oxford LABX®	32	28	36	22
SC-7332cap	Only cap, designed for Oxford LABX®	-	-	-	-
SC-8047	Cup, ring, single open ended, for Horiba® 47mm holder	42	36	19	18

One set contains 100 pieces. Caps and rings are packed separately! No cutting necessary!

Certified Reference Materials

Catalog	Description
	Petro standards liquid (ICP / XRF) This catalog contains: <ul style="list-style-type: none"> - sulfur in oil - sulfur in gasoline - lubricants - wear metals - metal concentrates - XRF oil package - XRF monitors
	Waterbased standards for AA and ICP This catalog gives an overview of liquid reference materials used for ICP / ICP-MS, IC and AA: <ul style="list-style-type: none"> - ICP / ICP-MS standards - EPA standards - Ion chromatography standards - Atomic absorption standards